

09853819 081001

H.A.T.S. - Patient Information

Reset: Caliper Color

Patient Name

Age Gender

Missing Teeth

Malformed Teeth

Doctor or Staff Performing Analysis

NOTE:
If a tooth is malformed or missing, it is recommended that a "zero" be entered into the analysis. The resulting correction necessary will indicate the size of the replacement tooth or teeth necessary to eliminate the discrepancy.

Figure 2A

H.A.T.S. - Patient Information

Reset: Caliper Color

Pat

Age Gender

Missing Teeth

Malformed Teeth

Doctor or Staff Performing Analysis

NOTE:
If a tooth is malformed or missing, it is recommended that a "zero" be entered into the analysis. The resulting correction necessary will indicate the size of the replacement tooth or teeth necessary to eliminate the discrepancy.

Figure 2B

H.A.T.S. Reset Calculate Print

Analysis Type
☒ 6 Tooth
☐ 12 Tooth

Input Tooth Sizes
☐ Input Sums

UPPER

TOOTH	3	4	5	6	7	8	9	10	11	12	13	14
AVERAGE SIZE	10.40	6.84	7.04	7.91	6.96	8.82	8.82	6.96	7.91	7.04	6.84	10.40
MEASURED SIZE	0	0	0	7.88	5.52	8.72	8.74	5.56	7.86	0	0	0

LOWER

MEASURED SIZE	0	0	0	6.86	5.87	5.36	5.38	5.88	6.84	0	0	0
AVERAGE SIZE	11.14	7.27	7.15	6.93	5.94	5.42	5.42	5.94	6.93	7.15	7.27	11.14
TOOTH	6	5	4	3	2	1	1	2	3	4	5	6
	30	29	28	27	26	25	24	23	22	21	20	19

Total Anterior Average Total Overall Average

Upper Sum 44.28 47.38 0 95.94

Lower Sum 36.19 36.58 0 87.70

<Back

Figure 2C

H.A.T.S. - Analysis Summary

Upper Sum 44.28 ~ 213

Lower Sum 36.19 ~ 214

Recommended Correction:

Increase upper tooth structure by 2.6 mm ~ 215

Reduce lower tooth structure by 2.01 mm ~ 216

G4C Patient Specific Arch Size SMALL ~ 217

Figure 2D




FIGURE 3




FIGURE 4